

PATARI Keyboard Shortcuts Guide

Using keyboard shortcuts can significantly speed up your analysis workflow. This is a short list containing some of the most useful shortcuts for traditional MSOT analysis, like navigating the viewer and annotating ROIs. A complete list of built-in napari shortcuts can be found and modified by clicking **napari** → **Preferences** → **Shortcuts**.

In the future, we can extend this list with custom shortcuts, depending on the analysis protocol.

Navigation & Layer Control

Browse through scan data, frames, and toggle between different image layers.

Action	Windows / Linux	macOS
Step through Frames / Wavelengths	Left / Right Arrow	Left / Right Arrow
Switch Active Slider (Frame vs. Wavelength)	Alt + Up / Down Arrow	Option + Up / Down Arrow
Switch Layer (Sequential)	Ctrl + Up / Down Arrow	Cmd + Up / Down Arrow
Switch Layer & Hide All Others (Solo View)	Shift + Alt + Up / Down Arrow	Shift + Option + Up / Down Arrow
Toggle Active Layer Visibility (Hide/Show)	V	V
Toggle Grid Mode	Ctrl + G	Cmd + G
Reset View to Original State	Ctrl + R	Cmd + R
Undo Last Action	Ctrl + Z	Cmd + Z
Quick Select Layer by Initial	First Letter of layer name (<i>Layer list must be selected</i>)	First Letter of layer name (<i>Layer list must be selected</i>)

ROI Annotation & Editing

Important: The ROIs layer must be actively selected in your layers list for these annotation shortcuts to function.

Action	Windows / Linux	macOS
Draw Ellipse ROI	E	E
Draw Rectangle ROI	R	R
Draw Polygon ROI	P	P
Finish Drawing Polygon	Double-Click	Double-Click
Modify Polygon Vertices	D	D
Pan & Zoom Mid-Draw (Hold key)	Spacebar	Spacebar
Select a ROI (Move/Resize mode)	S	S
Select All ROIs	A	A
Copy Selected ROI(s)	Ctrl + C	Cmd + C
Paste Selected ROI(s)	Ctrl + V	Cmd + V
Delete Selected ROI(s)	Delete OR Backspace	Delete OR Backspace
Save ROI Data*	Shift + Ctrl + S	Shift + Cmd + S

Image Processing

Action	Windows / Linux	macOS
Unmix using Default Preset*	Shift + Ctrl + U	Shift + Cmd + U

*not napari native